

Setup Guide for Brainify ML Flask App

0. Download the zip file from <http://detect-tumor.rf.gd/V1.zip>

1. Unzip the file

```
unzip V1.zip
```

2. Change the dir to Brainify

```
cd brainify
```

3. Install Dependencies

```
pip install -r requirements.txt
```

4. Setup Jupyter Notebook

```
cd notebooks
```

5. Run Jupyter Notebook

```
jupyter notebook
```

6. Prepare Training Data

Ensure that your training data images are located in the `data` folder within the notebook directory.

7. Training the Model

Follow the instructions in the notebook (`brainify_ml.ipynb`) to train and test your ML model using the provided data.

8. Deploying the Flask App

Navigate back to the root directory of your Brainify project and run the Flask app.

```
cd ...
```

```
python app.py
```

9. Accessing the Deployed Version <http://detect-tumor.rf.gd/>

Once the Flask app is running, the online version of your project can be accessed through the URL: detect-tumor.rf.gd.

10. Patient Login Credentials

To login as a patient on the website, use the following credentials:

- Email: patient@brainify.com
- Password: [123](#)